



Concentrated Quat Disinfectant

For Use In Hospitals, Food establishments, Industrial sites, and Institutions.

Description

This product is a one-step disinfectant, cleaner and deodorant designed for general cleaning, disinfecting, deodorizing and controlling mold and mildew on hard, non-porous non-food contact surfaces. This product is formulated for use on washable hard, non-porous surfaces such as: glass surfaces, aluminum, laminated surfaces, metal, glazed porcelain, glazed tile, sealed stone, plastic, chrome and vinyl, enameled surfaces, and painted woodwork.

Directions

FOR USE AS A ONE-STEP, GENERAL, HOSPITAL, MEDICAL DISINFECTANT, FUNGICIDE, VIRUCIDE*, CLEANER:

1. Pre-clean heavily soiled areas.
2. Apply use solution of 2 oz. of this product per gal. of water to disinfect hard, non-porous surfaces with a sponge, brush, cloth, mop, by immersion, auto scrubber, mechanical spray device, hand pump trigger spray device, coarse trigger spray device. For spray applications, spray 6-8 inches from surface. Do not breathe spray.
3. Treated surfaces must remain wet for 10 minutes.
4. Allow to air dry. Rinsing of floors is not necessary unless they are to be waxed or polished.
5. Prepare a fresh solution daily or when visibly dirty.

For detailed instructions and safety information, please review container label and SDS.



Features

- **Effective against a broad-spectrum of bacteria**
- **No-rinse neutral pH disinfectant**

Applications

One-step disinfectant for:

- Hospitals
- Food establishments
- Industrial sites
- Institutional kitchens

Specifications

ColorPurple
 Fragrance..... Floral
 pH Range 7.5 - 8.0
 Dilution 2 oz./gal.
 Density.....8.34
 EPA Registered Yes

Packaging

ITEM #	SIZE
BSL4740041	4/1 gallon



Contact your Account Representative for more information
 140 Private Brand Way, Athens, TN 37303 | mybrightsolutions.com
 In case of emergency call Infotrac 800.535.5053
 © 2018 Midlab, Inc. All right reserved.